

Sterilization Trays

Trays 20x10 and 20x15



ALUMINIUM TRAYS



Anodised aluminium - Colours



- Size : 20 x15 x 2 cm Ref. 181480
- Size : 20 x10 x 2 cm Ref. 181470




STAINLESS STEEL TRAYS



Stainless steel - High quality

Perfect for your exam sets.

- Size : 20 x15 x 2 cm Ref. 182480
- Size : 20 x10 x 2 cm Ref. 182470



ALUMINIUM RACKS

For 6 trays 20 x 15 :

- Size 16 x 20 x 18,5 cm Ref. 181485 (empty)

For 6 trays 20 x 10 :

- Size 10,3 x 20 x 18,5 cm Ref. 181475 (empty)



STAINLESS STEEL RACK

For 6 trays 20 x 15

- Size : 21 x 15,8 x 18,2 cm Ref. 182485 (empty)

For 6 trays 20 x 10

- Size : 20,6 x 10,2 x 18,2 cm Ref. 182475 (empty)

Can be wall mounted



Aluminium Tubs



Ultra light, for use as a portable operatory drawer. One Tub is recommended for each different type of procedure. Can be fitted with a divider and arranged in an aluminium rack.

- Aluminium Tub - Colours 1 2 3 4 - Size 33,5 x 23,5 x 5,2 cm Ref. 181510
- Tub Divider Ref. 181520
- Empty Rack (for 4 Tubs) - Size 34 x 24 x 41 cm Ref. 181525

Can be wall mounted

Sterilization Trays



Aluminium Steri Trays



Anodised aluminium - Colour **1** only. Multiperforated cover with lockable system. With replaceable paper filter. Silicone heat resistant lockable seal (230°). Aluminium instrument holder.



STERI TRAY 28X18

- Ref. 187450
- External size 29,5 x 19 x 4,6 cm
- Internal size 28 x 18 x 2,8 cm
- Ref. 187458 - Refill of 100 Paper Filters



STERI TRAY 18X14

- Ref. 187460
- External size 19 x 14,5 x 4,6 cm
- Internal size 18 x 13,5 x 2,8 cm
- Ref. 187468 - Refill of 100 Paper Filters



Sterilizing Containers

Anodised «high quality» aluminium - Colour **1** only. Multiperforated cover with lockable system. With replaceable paper filter. Silicone heat resistant seal (230°).



CONTAINER 25X15X6

- Ref. 300200
- External size 26 x 16 x 8,5 cm
- Internal size 25 x 15 x 5,8 cm
- Ref. 300250 - Refill of 100 Paper Filters



CONTAINER 25X10X5

- Ref. 300100
- External size 26 x 11 x 7,5 cm
- Internal size 25 x 10 x 5,1 cm
- Ref. 300150 - Refill of 100 Paper Filters



Steri Boxes

Highly efficient and ecological sterilization boxes! Fully airtight, they are equipped with a teflon filter effective up to 2000 cycles, a holder for an identification label and a seal system. Anodised aluminium - Colour **1** only. For suitable cassettes, please contact us.



STERI BOX 28X18

- Ref. 187440
- External size 30 x 19 x 6 cm
- Internal size 28 x 18 x 4 cm
- Ref. 187446 - Refill of Teflon Filter - per unit



Filtration system
with 2 silicone
gaskets and 1 Teflon
filter.



STERI BOX 18X14

- Ref. 187445
- External size 19 x 15,6 x 5,5 cm
- Internal size 18 x 13,5 x 3,5 cm
- Ref. 187447 - Refill of Teflon Filter - per unit

Sterilization Boxes

Aluminium Boxes



Sizes (cm)*	Grey Colour
17 x 7 x 3	180014-1
18 x 9 x 3	180016-1
18 x 9 x 4	180017-1
21 x 10 x 3	180019-1
21 x 10 x 4	180020-1
21 x 10 x 5	180021-1



Sizes (cm)*	Blue Colour	Red Colour	Green Colour	Yellow Colour
17 x 7 x 3	180014-2	180014-3	180014-4	180014-5
18 x 9 x 3	180016-2	180016-3	180016-4	180016-5
21 x 10 x 3	180019-2	180019-3	180019-4	180019-5
21 x 10 x 5	180021-2	180021-3	180021-4	180021-5

Aluminium Boxes Small Sizes



BOXES 5 X 4 X 3 CM*

Colour	Reference
Grey Colour	180004-1
Blue Colour	180004-2
Red Colour	180004-3
Green Colour	180004-4
Yellow Colour	180004-5

BOXES 10 X 3 X 2 CM*

Colour	Reference
Grey Colour	180005-1
Blue Colour	180005-2
Red Colour	180005-3
Green Colour	180005-4

BOXES 13 X 5 X 3 CM*

Colour	Reference
Grey Colour	180007-1
Blue Colour	180007-2
Red Colour	180007-3
Green Colour	180007-4


* Above mentioned sizes correspond to our trade descriptions. Precise dimensions are available on request.

Sterilization Boxes

Alu Perforated Boxes



Colour **1** only



SPECIAL AUTOCLAVE

Sizes (cm)*	Reference
17 x 7 x 3	184014
18 x 9 x 3	184016
21 x 10 x 3	184019
21 x 10 x 5	184021

Stainless Steel Boxes



Sizes (cm)*	Reference
17 x 7 x 3	182014
18 x 9 x 3	182016
21 x 10 x 3	182019
21 x 10 x 5	182021

Stainless Steel Perforated Boxes



SPECIAL AUTOCLAVE

Sizes (cm)*	Reference
17 x 7 x 3	185014
18 x 9 x 3	185016
21 x 10 x 3	185019
21 x 10 x 5	185021

Alu Boxes Large Sizes

High quality anodised aluminium, 15/10 thickness.

Colour **1** only.   



Sizes (cm)*	Reference
25 x 10 x 5	180022-1
25 x 15 x 6	180024-1

Polybox




Colour	Not Perforated	Perforated
Blue 2	410016-2	412016-2
Green 4	410016-4	412016-4
Ivory 1	410016-1	412016-1

Size 18 x 9 x 3 cm* - Polycarbonate, very light and resistant material. **Transparent cover allows direct view of the instruments.**

* Above mentioned sizes correspond to our trade descriptions. Precise dimensions are available on request.

Retractors



LOWER THIRD MOLAR RETRACTOR

- Ref. 074000 - The pair (right + left)

Specifically designed for removal of lower wisdom teeth. Claws for support. Improve visibility and access to the third molar region, by repelling the labial commissure both horizontally and vertically. Very comfortable handling.



UPPER THIRD MOLAR RETRACTOR

- Ref. 071700 - The pair (right + left)

Specifically designed for removal of upper wisdom teeth. Claws for support. Improve visibility and access to the third molar region. The important factor is that the flap is pushed instead of being pulled.



FLAP RETRACTOR

- Ref. 071108 - Per unit, width 8 mm
- Ref. 071112 - Per unit, width 12 mm

For removal of lower or upper wisdom teeth. Claws for support - Can be used for removing the upper or lower wisdom teeth, also for all types of oral surgery.



GINESTET TYPE RETRACTOR

- Ref. 072400 - The pair (right + left)

Keep the flap retracted during surgeries on impacted wisdom teeth.



DAUTREY TYPE RETRACTOR

- Ref. 072200 - Per unit, short (H. 6 et 4,5 cm)
- Ref. 072201 - Per unit, long (H. 9 et 7,5 cm)

Those retractors are curved, with very smooth edges, so they don't hurt commissures during surgeries on wisdom teeth.



FARABEUZ TYPE RETRACTOR

- Ref. 072100 - Set of 2

This retractor is multi-purpose, it can be used for all types of oral surgery.

Retractors and Mouth Props



MOUTH PROPS

- Ref. 070300 - Set of 2 (in an aluminium box)
- 2 sizes : 15 and 25 mm

Non reflective stainless steel. The patient bites on a soft silicone surface to keep the mouth open. Available in 2 sizes : 15 and 25 mm.



«VISIO-N» RETRACTOR

- Ref. 075100 - Per unit

Convenient to all types of daily procedures, creates the perfect gap between the cheek and jaw, comfortable for the patient as well as for the user.



MALLEABLE RETRACTORS

- Ref. 071000 - Set of 3
- 3 widths : 10, 15 and 20 mm


Can be easily bent with fingers or pliers to the shape you require.



UNIVERSAL RETRACTOR

- Ref. 075000 - Per unit

For all types of oral surgery.



RETRACTOR FOR LINGUAL APPROACH (1)

- Ref. 070800 - Set of 2

This retractor is designed for the lingual approach of the lower wisdom teeth.

TONGUE RESTRAINER (2)

- Ref. 072000 - Set of 2

This retractor allows you to gently keep the tongue down.




"S" SHAPE RETRACTOR

- Ref. 072500 - Per unit

Curved edges, retract both a flap and cheek without trauma. 2 ends of different widths : 12 and 20 mm.

STAINLESS STEEL CHEEK RETRACTORS




REVERSIBLE RETRACTOR

- Ref. 072600 - Set of 2
- One large + one small end.



SELF-RETAINING RETRACTORS

- Ref. 072753 - Small
- Ref. 072754 - Large
- Keep both cheeks and lips retracted without trauma.



Various Accessories



TONGUE PROTECTOR

- Ref. 190514 - The pair (right + left)


Without mirror. Tongue and cheek protection.
Easy access to the whole arch. Allows you to use a handpiece at any angle safely.



MIRROR TONGUE PROTECTOR

- Ref. 190515 - Per unit


With a high quality mirror. Tongue and cheek protection.
Easy access to the whole arch. Allows you to use a handpiece at any angle safely.



EXAMINATION MIRRORS

- Ref. 193510 - Set of 6

Anodised aluminium handles - Mirror Ø 2 cm. Sold as set of 6 with assorted coloured handles : 1 2 3 4 5 6.



PARALLELISM MIRROR

- Ref. 190750 - Sold with handle

High Quality front surface mirror with laser marked lines on the surface. **Parallel lines are spaced 3 mm, making it easy to check abutments parallelism.**
Size of the mirror 4,5 x 2,9 cm



MIRRORS Ø4 & Ø5 CM

- Ref. 190772 - Mirror Ø 4 cm
- Ref. 190773 - Mirror Ø 5 cm

For Intra-Oral Photography or for Motivation.
Sold with handle.



SILICONE CODE RING

- Ref. 190202 - Box of 50 pieces

4 colours 2 3 4 5 - For identifying and organising instruments. **Autoclavable 134°C.**

Various Accessories



BOWLS

- **Small Size** - Ref. 182041
Ø 5,3 cm - H. 2,5 cm - capacity 1,8 cl
 - **Large Size** - Ref. 182042
Ø 6,1 cm - H. 2,7 cm - capacity 3,4 cl
- Stainless Steel.
Small bowls for rinse, dilute and mix.



COTTON HOLDERS

- Size Ø 5,5 cm - H. 7,5 cm
- Ref. 190300-1 - Alu colour 1
 - Ref. 190300-2 - Alu colour 2
 - Ref. 190300-4 - Alu, colour 4
 - Ref. 192300 - Stainless Steel



GERMICIDE BOXES

- Ref. 182050 - Size Ø 7,5 cm - H. 3,5 cm
- Stainless Steel. Composed of 3 parts : 1 bowl, 1 perforated insert and 1 opening cover with slot.
Prevent airborne contamination.



WASTE BINS

- Size Ø 7 cm - H. 8 cm
- Ref. 190400-1 Alu colour 1
 - Ref. 190400-2 Alu colour 2
 - Ref. 190400-4 Alu colour 4
 - Ref. 192400 - Stainless Steel
- Inside bowl



BITE REGISTRATION TRAYS RIGID

- Ref. 190150 - Set of 2
- Stainless steel. This bite is designed to compress the impression material.



WIRE BITE TRAYS

- Ref. 190110 - Box of 10
- Stainless steel. This wire bite makes the handling of the impression material easier.

Silicone Refills

PLUGS FOR BUR STANDS

- Ref. 190200
- Set of 8



Sold as a set of 8.

Colours 1 2 3 4 5 6 7 8 9 10 11 12 13.

Used for Plug In, Ultralight Bur Holders (p.10-11-12), Endo Star (p.15) and Implantology Kits (p.22-23).

PLUGS FOR ENDODONTIC STANDS

- Ref. 182175
- Set of 8



Sold as a set of 8.


Colour 2 only.

Used for Endo Star and Endo Ultralight (p.15).

SILICONES TUBES

Set of 6, colours 2 3 4 5

Central slot is available in 3 different sizes :



- Ref. 500600 - Set of 6 tubes
- For instruments Ø 2 à 3,5 mm



- Ref. 500610 - Set of 6 tubes
- For instruments Ø 3,5 à 5 mm



- Ref. 500620 - Set of 6 tubes
- For instruments Ø 5 à 6,5 mm

For Implantology Kits (p.22-23), Multi Holders, Multi Clip Holders, and Multi Plug In (p.12).

CURSORS FOR STERIMETER

- Ref. 182370
- Set of 12




Colours 2 3 4 5.

For all endodontic stands (except Aluminium range) (p.14-15)

INDICATOR STRIPS

- Ref. 190275
- Set of 12



Colours 2 3 4 5.

For Bur Clips (p.8), Multi Plugs (p.12), Bur Flash (p.13), Micro Cassettes and Divided Cassettes (p.6).

FOR ERGO CLIPS / ERGO CASSETTES (page 3) and FLEXI CLIPS (off catalogue)

Colours : 1 2 3 4 5 6 7 8 9 13



Size and Cassette Ref.		Composition	Refill Ref.
18x8 cm	182703/182803	3 notches ø 10 mm	182713
18x11 cm	182700/182753/182800	5 notches ø 8,5 mm	182705
	182801	5 notches ø 11 mm	182724
18x15 cm	182710/182750/182800	10 notches ø 8,5 mm	182715
	182811	10 notches ø 11 mm	182729
23x18 cm	182712/182752	10 notches ø 8,5 mm + flat strip	182722
	182812	15 notches ø 8,5 mm	182743
28x18 cm	182720/182751/182820	20 notches ø 8,5 cm	182736
	182725	10 notches ø 8,5 mm + waves	182721



FOR FLAT CASSETTES (page 3) - Colours : 1 2 3 4 5 6 7 8 9 13



Size and Cassette Ref.		Composition	Refill Ref.
18x6 cm	182063/183063	5 notches ø 8,5 mm	182073
18x9 cm	182055/183055	7 notches ø 8,5 mm	182089
18x13 cm	182062/183062	10 notches ø 8,5 mm	182072



Silicone Refills


FOR CLASSIC STYLE & GALAXY CASSETTES (page 2) AND ULTRALIGHT CASSETTES (off catalogue)

Colours : 2 3 4 5



Size and Cassette Ref.		Composition	Refill Ref.
Size 18x7	182069	4 notches ø 11 mm	182078
	182039	5 notches ø 8 mm	182052
	182079	2 hand pieces notches	182088
	182044	5 notches Galaxy	182045
Size 18x10	182054	6 notches ø 11 mm	182084
	182025	7 notches Galaxy	182026
Size 18x14	182066	8 notches ø 11 mm	182077
	182056/182036	10 notches ø 8 cm	182083
	182046	10 notches Galaxy	182047
Size 18x20	182057	8 notches ø 11 mm + clear area	182082
	182048	10 notches Galaxy + clear area	182028
Size 28x18	182068	16 notches ø 11 mm	182071
	182058	20 notches ø 8 cm	182081
	182049	20 notches Galaxy	182029
	182087	10 notches ø 11 mm + clear area	182074
	182037	12 notches ø 8 mm + clear area	182059
	182027	12 notches Galaxy + clear area	182024

FOR EASY CLIP CASSETTES (page 4) - Colours : 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.		Composition	Refill Ref.
18x5 cm	182903	3 notches ø 10 mm	182953
	182904	4 notches ø 9 mm	182954
18x8 cm	182900	5 notches ø 8,5 mm	182955
	182901	5 notches ø 11 mm	182956
18x13 cm	182910	10 notches ø 8,5 mm	182960
	182911	10 notches ø 11 mm	182961
20x18 cm	182912	15 notches ø 8,5 mm	182965
25x18 cm	182920	20 notches ø 8,5 cm	182970

FOR EASY TRAY CASSETTES (page 5) - Colours 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.		Composition	Refill Ref.
18x5 cm	183703	3 notches ø 10 mm	182953
	183704	4 notches ø 9 mm	182954
18x7 cm	183700	5 notches ø 8,5 mm	183701
18x10 cm	183705	5 notches ø 8,5 mm + flat strip	183706
18x14 cm	183710	10 notches ø 8,5 mm + flat strip	183711
21x18 cm	183715	15 notches ø 8,5 cm + flat strip	183716
28x18 cm	183720	20 notches ø 8,5 cm + flat strip	183721

FOR EASY TRAY ORTHO CASSETTES (page 7) - Colours 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.		Composition	Refill Ref.
18x14 cm	183730	2 notches ø 8,5 mm + waves	183731
21x18 cm	183732	3 notches ø 8,5 mm + waves	183733
28x18 cm	183734	3 notches ø 8,5 mm + waves	183735

La créativité de l'artisan
The expertise of craftsmen

La rigueur de l'industriel
Industrial Excellence



Distributeur / Dealer :

Pictogrammes : Bain à Ultrasons  Thermo-laveur  Autoclave  134°C 

NICHROMINOX - 18 rue des Frères Lumière - F-69720 SAINT-BONNET-DE-MURE - FRANCE
 Tel. +33 (0)4 78 74 04 15 - Fax +33 (0)4 78 01 61 80 - email : office@nichrominox.fr - www.nichrominox.fr

Valid from 01.01.2019 to 31.12.2019